


LeviKam 8 point Sling

L8-U-1



The LeviKam eight-point transfer sling is compatible with hoists for people with disabilities. It features high-quality workmanship and a well-thought-out design. It is made of durable, safe polyester mesh that is highly water-permeable and dries quickly. The sling allows users to be transported in a lying position, ensuring stability and comfort. The user's position can be adjusted to their needs thanks to the eight-point, adjustable suspension.



*The values given are approximate due to differences in body structure.

Main advantages:

- Maximum load capacity: 200 kg,
- Eight attachment points enable comfortable and stable transfers in a lying position,
- Multi-level adjustment of the sling attachments,
- The slings are equipped with a physiological cutout that facilitates toileting and hygiene procedures,
- The seats are available in two versions: a standard bath seat (made of polyester mesh) or a special comfort seat (made of solid polyester material, sewn upon request).

Instructional video:

Putting on vests in a lying position:



Using the 8-point sling for transfer:

