

LeviKam Active with head support Sling

LACPL-XS/S/M/L/XL-1



The LeviKam Active sling with head support is compatible with hoists for people with disabilities. It features high-quality workmanship and a well-thought-out design. It is made of durable, safe polyester mesh that is highly water-permeable and dries quickly. The sling allows users to transport the user in a seated position while providing back support. The user's position can be adjusted to their needs thanks to the four-point, adjustable suspension. A sewn-in maneuvering handle allows the user to be rotated while hanging.



LeviKam Active with Head Support					
	XS size	S size	M size	L size	XL size
A size	84	86	96	107	122
B size	92	94	96	98	115
C size	38	40	41	44	45

*The values given are approximate due to differences in body structure.

Main advantages:

- Maximum load capacity: up to 200 kg,
- Wide range of sizes available: from XS to XL,
- Multi-level adjustment of the sling hooks,
- Slings feature a physiological cutout to facilitate toileting and hygiene procedures,
- Seats are available in two versions: a standard bath seat (made of polyester mesh) or a special comfort seat (made of solid polyester material, sewn upon request).

Instructional videos:

Putting on sling in a sitting position:



Putting on sling in a lying position:

