

# LeviKam Plus LPL-XS/S/M/L/XL-1

## Sling



The LeviKam Plus transfer sling is compatible with hoists for people with disabilities. It features high-quality workmanship and a well-thought-out design. It is made of durable, safe polyester mesh that is highly water-permeable and dries quickly. The LeviKam Plus sling allows the user to be transported in a seated position while providing back support. The user's position can be adjusted to their needs thanks to the six-point, adjustable suspension. A sewn-in maneuvering handle allows the user to be rotated while hanging.



LeviKam +	XS size	S size	M size	L size	XL size
A size	66	79	90	101	114
B size	34	37	44	46	48
C size	50	53	52	54	57
D size	51	56	58	75	88

\*The values given are approximate due to differences in body structure.

### Main advantages of slings:

- Maximum load capacity: 200 kg,
- Head support,
- Multi-level adjustment of the sling attachment points,
- Wide range of sizes available: from XS to XL,
- Slings feature a physiological cutout to facilitate toileting and hygiene procedures,
- Seats are available in two versions: a standard bath seat (made of polyester mesh) or a special comfort seat (made of solid polyester material, sewn upon request).



### Instructional videos

Putting on slings in a sitting position:



Putting on slings in a lying position:

